

CAB22004

Product Information

| | | |
|------------------------------|-----------------------------|------------------------------------|
| Product SKU: CAB22004 | Gene ID: 9589 | Size: 20uL, 100uL |
| Clone No: - | Host Species: Rabbit | Reactivity: Human,Mouse,Rat |

Additional Information

| | |
|-----------------------------|--------------------------------|
| Observed MW: 55kDa | Conjugate: Unconjugated |
| Calculated MW: 44kDa | Isotype: IgG |

Immunogen Information

Background: The Wilms tumor suppressor gene WT1 appears to play a role in both transcriptional and posttranscriptional regulation of certain cellular genes. This gene encodes a WT1-associating protein, which is a ubiquitously expressed nuclear protein. Like WT1 protein, this protein is localized throughout the nucleoplasm as well as in speckles and partially colocalizes with splicing factors. Alternative splicing of this gene results in several transcript variants encoding three different isoforms.

Recommended Dilution: WB,1:100 - 1:500 IF/ICC,1:50 - 1:200

Synonyms: Mum2; WTAP

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 212-355 of human WTAP (NP_004897.2).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.