

CAB22017

Product Information

Product SKU:	CAB22017	Gene ID:	6423	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	33kDa	Conjugate:	Unconjugated
Calculated MW:	33kDa	Isotype:	IgG

Immunogen Information

Background: This gene encodes a member of the SFRP family that contains a cysteine-rich domain homologous to the putative Wnt-binding site of Frizzled proteins. SFRPs act as soluble modulators of Wnt signaling. Methylation of this gene is a potential marker for the presence of colorectal cancer.

Recommended Dilution: WB,1:500 - 1:1000 IF/ICC,1:50 - 1:200

Synonyms: FRP-2; SARP1; SDF-5; [KD Validated] SFRP2

Purification Method: Affinity purification

Immunogen: A synthetic peptide corresponding to a sequence within amino acids 1-100 of human SFRP2 (NP_003004.1).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.