

## CAB22022

---

### Product Information

<b>Product SKU:</b> CAB22022	<b>Gene ID:</b> 3308	<b>Size:</b> 20uL, 100uL
<b>Clone No:</b> -	<b>Host Species:</b> Rabbit	<b>Reactivity:</b> Human,Mouse,Rat

---

### Additional Information

<b>Observed MW:</b> 110kDa	<b>Conjugate:</b> Unconjugated
<b>Calculated MW:</b> 94kDa	<b>Isotype:</b> IgG

---

### Immunogen Information

**Background:** Predicted to enable ATP binding activity. Involved in chaperone-mediated protein complex assembly and protein insertion into mitochondrial outer membrane. Located in cytosol and extracellular exosome. Implicated in Chagas disease. Biomarker of chronic obstructive pulmonary disease; rheumatoid arthritis; type 2 diabetes mellitus; and ulcerative colitis.

**Recommended Dilution:** WB,1:500 - 1:1000

**Synonyms:** RY; APG-2; HSPH2; hsp70; hsp70RY; HEL-S-5a; HS24/P52; HSPA4

**Purification Method:** Affinity purification

**Immunogen:** Recombinant fusion protein containing a sequence corresponding to amino acids 680-840 of human HSPA4 (NP\_002145.3).

**Storage:** Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.