

## CAB22036

---

### Product Information

<b>Product SKU:</b>	CAB22036	<b>Gene ID:</b>	5155	<b>Size:</b>	20uL, 100uL
<b>Clone No:</b>	ARC53713	<b>Host Species:</b>	Rabbit	<b>Reactivity:</b>	Mouse

---

### Additional Information

<b>Observed MW:</b>	Refer to figures	<b>Conjugate:</b>	Unconjugated
<b>Calculated MW:</b>	27kDa	<b>Isotype:</b>	IgG

---

### Immunogen Information

<b>Background:</b>	This gene encodes a member of the protein family comprised of both platelet-derived growth factors (PDGF) and vascular endothelial growth factors (VEGF). The encoded preproprotein is proteolytically processed to generate platelet-derived growth factor subunit B, which can homodimerize, or alternatively, heterodimerize with the related platelet-derived growth factor subunit A. These proteins bind and activate PDGF receptor tyrosine kinases, which play a role in a wide range of developmental processes. Mutations in this gene are associated with meningioma. Reciprocal translocations between chromosomes 22 and 17, at sites where this gene and that for collagen type 1, alpha 1 are located, are associated with dermatofibrosarcoma protuberans, a rare skin tumor. Alternative splicing results in multiple transcript variants.
<b>Recommended Dilution:</b>	IF/ICC,1:50 - 1:200
<b>Synonyms:</b>	SIS; SSV; IBGC5; PDGF2; c-sis; PDGF-2; PDGFB
<b>Purification Method:</b>	Affinity purification
<b>Immunogen:</b>	A synthetic peptide corresponding to a sequence within amino acids 142-241 of human PDGFB (NP_002599.1).
<b>Storage:</b>	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.05% proclin300,0.05% BSA,50% glycerol,pH7.3.

**Contact Details | Dublin, Ireland**

**Email:** [techsupport@assaygenie.com](mailto:techsupport@assaygenie.com) | **Web:** [www.assaygenie.com](http://www.assaygenie.com)

Copyright © 2020 Reagent Genie, All Rights Reserved. All information / detail is correct at time of going to print.