

CAB22049

Product Information

Product SKU:	CAB22049	Gene ID:	9340	Size:	20uL, 100uL
Clone No:	ARC54969	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	70kDa	Conjugate:	Unconjugated
Calculated MW:	63kDa	Isotype:	IgG

Immunogen Information

Background: This gene encodes a G protein-coupled receptor that is closely related to the glucagon receptor and binds to glucagon-like peptide-2 (GLP2). Signalling through GLP2 stimulates intestinal growth and increases villus height in the small intestine, concomitant with increased crypt cell proliferation and decreased enterocyte apoptosis.

Recommended Dilution: WB,1:2000 - 1:20000

Synonyms: GLP2R

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 460-553 of human GLP2R (NP_004237.1).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.05% proclin300,0.05% BSA,50% glycerol,pH7.3.