

CAB22057

Product Information

Product SKU: CAB22057	Gene ID: 55177	Size: 20uL, 100uL
Clone No: ARC54933	Host Species: Rabbit	Reactivity: Rat

Additional Information

Observed MW: 52kDa	Conjugate: Unconjugated
Calculated MW: 52kDa	Isotype: IgG

Immunogen Information

Background: Enables microtubule binding activity. Involved in cellular calcium ion homeostasis. Located in several cellular components, including intercellular bridge; mitochondrial outer membrane; and spindle.

Recommended Dilution: WB, 1:1000 - 1:5000

Synonyms: RMD3; RMD-3; FAM82C; FAM82A2; ptpip51; PTPIP51/RMDN3

Purification Method: Affinity purification

Immunogen: A synthetic peptide corresponding to a sequence within amino acids 200-300 of human PTPIP51/RMDN3 (NP_060615.1).

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300, 0.05% BSA, 50% glycerol, pH 7.3.