

## CAB22106

---

### Product Information

<b>Product SKU:</b>	CAB22106	<b>Gene ID:</b>	64326	<b>Size:</b>	20uL, 100uL
<b>Clone No:</b>	-	<b>Host Species:</b>	Rabbit	<b>Reactivity:</b>	Human

---

### Additional Information

<b>Observed MW:</b>	102kDa	<b>Conjugate:</b>	Unconjugated
<b>Calculated MW:</b>	80kDa	<b>Isotype:</b>	IgG

---

### Immunogen Information

**Background:** Enables ubiquitin protein ligase activity. Involved in positive regulation of proteasomal ubiquitin-dependent protein catabolic process; proteasome-mediated ubiquitin-dependent protein catabolic process; and response to ionizing radiation. Part of Cul4A-RING E3 ubiquitin ligase complex.

**Recommended Dilution:** WB, 1:500 - 1:1000

**Synonyms:** FAP78; RFWD2; CFAP78; RNF200; COP1

**Purification Method:** Affinity purification

**Immunogen:** Recombinant fusion protein containing a sequence corresponding to amino acids 582-731 of human COP1 (NP\_071902.2).

**Storage:** Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300, 50% glycerol, pH7.3.