

CAB22114

Product Information

Product SKU:	CAB22114	Gene ID:	12525	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Mouse

Additional Information

Observed MW:	35-42kDa	Conjugate:	Unconjugated
Calculated MW:	27kDa	Isotype:	IgG

Immunogen Information

Background: Enables identical protein binding activity. Acts upstream of or within several processes, including cytotoxic T cell differentiation; defense response to virus; and positive regulation of calcium-mediated signaling. Located in external side of plasma membrane. Is expressed in several structures, including endocrine gland; exocrine gland; genitourinary system; gut; and hemolymphoid system gland. Orthologous to human CD8A (CD8a molecule).

Recommended Dilution: WB, 1:500 - 1:1000

Synonyms: Ly-2; Ly-B; Ly-35; Lyt-2; CD8A

Purification Method: Affinity purification

Immunogen: Recombinant protein of mouse CD8A

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300, 50% glycerol, pH7.3.