

CAB22116

Product Information

Product SKU: CAB22116	Gene ID: 84879	Size: 20uL, 100uL
Clone No: -	Host Species: Rabbit	Reactivity: Human,Mouse,Rat

Additional Information

Observed MW: 73kDa	Conjugate: Unconjugated
Calculated MW: 60kDa	Isotype: IgG

Immunogen Information

Background: The protein encoded by this gene is a transmembrane protein and sodium-dependent lysophosphatidylcholine transporter. The encoded protein is involved in the establishment of the blood-brain barrier and is required for brain growth and function. Defects in this gene are a cause of a progressive microcephaly syndrome.

Recommended Dilution: WB,1:500 - 1:1000

Synonyms: NLS1; MFSD2; MCPH15; SLC59A1; HsMFSD2A; NEDMISBA; MFSD2A

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence within amino acids 464-543 of human MFSD2A (NP_001129965.1).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.