

CAB22171

Product Information

Product SKU:	CAB22171	Gene ID:	969	Size:	20uL, 100uL
Clone No:	ARC56158	Host Species:	Rabbit	Reactivity:	Mouse,Rat

Additional Information

Observed MW:	25-45kDa	Conjugate:	Unconjugated
Calculated MW:	23kDa	Isotype:	IgG

Immunogen Information

Background: This gene encodes a member of the calcium dependent lectin superfamily of type II transmembrane receptors. Expression of the encoded protein is induced upon activation of T lymphocytes, and may play a role in proliferation. Furthermore, the protein may act to transmit signals in natural killer cells and platelets.

Recommended Dilution: WB,1:2000 - 1:4000 IF/ICC,1:50 - 1:200

Synonyms: AIM; EA1; MLR-3; CLEC2C; GP32/28; BL-AC/P26; CD69

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 62-199 of human CD69 (NP_001772.1).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.05% proclin300,0.05% BSA,50% glycerol,pH7.3.