

## CAB2219

---

### Product Information

<b>Product SKU:</b>	CAB2219	<b>Gene ID:</b>	8863	<b>Size:</b>	20uL, 100uL
<b>Clone No:</b>	ARC1893	<b>Host Species:</b>	Rabbit	<b>Reactivity:</b>	Human,Mouse,Rat

---

### Additional Information

<b>Observed MW:</b>	132kDa	<b>Conjugate:</b>	Unconjugated
<b>Calculated MW:</b>	132kDa	<b>Isotype:</b>	IgG

---

### Immunogen Information

**Background:** This gene is a member of the Period family of genes and is expressed in a circadian pattern in the suprachiasmatic nucleus, the primary circadian pacemaker in the mammalian brain. Genes in this family encode components of the circadian rhythms of locomotor activity, metabolism, and behavior. This gene is upregulated by CLOCK/ARNTL heterodimers but then represses this upregulation in a feedback loop using PER/CRY heterodimers to interact with CLOCK/ARNTL. Polymorphisms in this gene have been linked to sleep disorders. Multiple transcript variants encoding different isoforms have been found for this gene.

**Recommended Dilution:** WB,1:500 - 1:1000

**Synonyms:** GIG13; FASPS3; PER3

**Purification Method:** Affinity purification

**Immunogen:** Recombinant fusion protein containing a sequence corresponding to amino acids 10-100 of human PER3 (P56645).

**Storage:** Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,0.05% BSA,50% glycerol,pH7.3.