

CAB2222

Product Information

Product SKU:	CAB2222	Gene ID:	51692	Size:	20uL, 100uL
Clone No:	ARC2567	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	75kDa	Conjugate:	Unmodified
Calculated MW:	77kDa	Isotype:	IgG

Immunogen Information

Background: This gene encodes a member of the metallo-beta-lactamase family. The encoded protein is a 73kDa subunit of the cleavage and polyadenylation specificity factor and functions as an endonuclease that recognizes the pre-mRNA 3'-cleavage site AAUAAA prior to polyadenylation. It also cleaves after the pre-mRNA sequence ACCCA during histone 3'-end pre-mRNA processing.

Recommended Dilution: WB,1:500 - 1:1000

Synonyms: CPSF73; NEDMHS; CPSF-73; NEDMHSN

Purification Method: Affinity purification

Immunogen: A synthetic peptide corresponding to a sequence within amino acids 585-684 of human CPSF3 (Q9UKF6).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,0.05% BSA,50% glycerol,pH7.3.