

CAB22233

Product Information

Product SKU: CAB22233	Gene ID: 2804	Size: 20uL, 100uL
Clone No: -	Host Species: Rabbit	Reactivity: Human

Additional Information

Observed MW: 376kDa	Conjugate: Unconjugated
Calculated MW: 376kDa	Isotype: IgG

Immunogen Information

Background: Enables RNA binding activity. Involved in protein localization to pericentriolar material. Located in Golgi apparatus and endoplasmic reticulum-Golgi intermediate compartment.

Recommended Dilution: WB, 1:500 - 1:1000 IF/ICC, 1:50 - 1:200

Synonyms: GCP; GCP372; GOLIM1; Giantin

Purification Method: Affinity purification

Immunogen: Recombinant protein of human Giantin

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300, 50% glycerol, pH7.3.