

CAB22236

Product Information

Product SKU: CAB22236	Gene ID: 11100	Size: 20uL, 100uL
Clone No: -	Host Species: Rabbit	Reactivity: Human, Mouse

Additional Information

Observed MW: 96kDa	Conjugate: Unconjugated
Calculated MW: 83kDa	Isotype: IgG

Immunogen Information

Background: This gene encodes a nuclear RNA-binding protein of the heterogeneous nuclear ribonucleoprotein (hnRNP) family. This protein binds specifically to adenovirus early-1B-55kDa oncoprotein. It may play an important role in nucleocytoplasmic RNA transport, and its function is modulated by early-1B-55kDa in adenovirus-infected cells.

Recommended Dilution: WB, 1:100 - 1:500

Synonyms: E1BAP5; E1B-AP5; HNRNPUL1; HNRNPUL1

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 209-410 of human HNRNPUL1 (NP_008971.2).

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300, 50% glycerol, pH7.3.