

## CAB22237

---

### Product Information

<b>Product SKU:</b> CAB22237	<b>Gene ID:</b> 11319	<b>Size:</b> 20uL, 100uL
<b>Clone No:</b> -	<b>Host Species:</b> Rabbit	<b>Reactivity:</b> Human,Mouse,Rat

---

### Additional Information

<b>Observed MW:</b> 76kDa	<b>Conjugate:</b> Unconjugated
<b>Calculated MW:</b> 73kDa	<b>Isotype:</b> IgG

---

### Immunogen Information

**Background:** Enables histone acetyltransferase binding activity. Involved in positive regulation of transcription by RNA polymerase II. Located in cytosol and nucleoplasm.

**Recommended Dilution:** WB,1:500 - 1:1000

**Synonyms:** GCR2; SGT1; HSGT1; ECD

**Purification Method:** Affinity purification

**Immunogen:** Recombinant fusion protein containing a sequence corresponding to amino acids 435-644 of human ECD (NP\_009196.1).

**Storage:** Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.