

CAB22251

Product Information

Product SKU:	CAB22251	Gene ID:	829547	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Arabidopsis thaliana

Additional Information

Observed MW:	50kDa	Conjugate:	Unconjugated
Calculated MW:	-	Isotype:	IgG

Immunogen Information

Background: Encodes an ABA-responsive element-binding protein with similarity to transcription factors that is expressed in response to stress and abscisic acid.

Recommended Dilution: WB, 1:100 - 1:500

Synonyms: abscisic acid responsive elements-binding factor 3; AtABF3; DC3 PROMOTER-BINDING FACTOR 5; DPBF5; F17I5.200; ABF3

Purification Method: Affinity purification

Immunogen: Recombinant protein of Arabidopsis thaliana ABF3

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300, 50% glycerol, pH7.3.