

CAB22258

Product Information

Product SKU:	CAB22258	Gene ID:	388	Size:	20uL, 100uL
Clone No:	ARC56483	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	22kDa	Conjugate:	Unconjugated
Calculated MW:	22kDa	Isotype:	IgG

Immunogen Information

Background: Predicted to enable GTP binding activity; GTPase activity; and protein kinase binding activity. Involved in several processes, including cellular response to hydrogen peroxide; cellular response to ionizing radiation; and regulation of cell migration. Located in cleavage furrow and endosome membrane. Biomarker of breast cancer.

Recommended Dilution: WB,1:2000 - 1:20000

Synonyms: ARH6; ARHB; RHOH6; MST081; MSTP081; RhoB

Purification Method: Affinity purification

Immunogen: A synthetic peptide corresponding to a sequence within amino acids 70-150 of human RhoB (NP_004031.1).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.05% proclin300,0.05% BSA,50% glycerol,pH7.3.