

CAB22275

Product Information

Product SKU:	CAB22275	Gene ID:	7465	Size:	20uL, 100uL
Clone No:	ARC57004	Host Species:	Rabbit	Reactivity:	Mouse

Additional Information

Observed MW:	Refer to figures	Conjugate:	Unconjugated
Calculated MW:	72kDa	Isotype:	IgG

Immunogen Information

Background: This gene encodes a nuclear protein, which is a tyrosine kinase belonging to the Ser/Thr family of protein kinases. This protein catalyzes the inhibitory tyrosine phosphorylation of CDC2/cyclin B kinase, and appears to coordinate the transition between DNA replication and mitosis by protecting the nucleus from cytoplasmically activated CDC2 kinase.

Recommended Dilution: IHC-P,1:50 - 1:200

Synonyms: WEE1A; WEE1hu; WEE1

Purification Method: Affinity purification

Immunogen: A synthetic peptide corresponding to a sequence within amino acids 1-100 of human WEE1 (NP_003381.1).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.05% proclin300,0.05% BSA,50% glycerol,pH7.3.