

CAB22276

Product Information

Product SKU:	CAB22276	Gene ID:	965	Size:	20uL, 100uL
Clone No:	ARC54834	Host Species:	Rabbit	Reactivity:	Human

Additional Information

Observed MW:	55-70kDa	Conjugate:	Unconjugated
Calculated MW:	28kDa	Isotype:	IgG

Immunogen Information

Background: This gene encodes a member of the immunoglobulin superfamily. The encoded protein is a ligand of the T lymphocyte CD2 protein, and functions in adhesion and activation of T lymphocytes. The protein is localized to the plasma membrane. Alternatively spliced transcript variants have been described.

Recommended Dilution: WB, 1:2000 - 1:20000

Synonyms: ag3; LFA3; LFA-3; CD58

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 29-215 of human CD58 (NP_001770.1).

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300, 0.05% BSA, 50% glycerol, pH 7.3.