

## CAB22286

---

### Product Information

<b>Product SKU:</b>	CAB22286	<b>Gene ID:</b>	3556	<b>Size:</b>	20uL, 100uL
<b>Clone No:</b>	ARC56983	<b>Host Species:</b>	Rabbit	<b>Reactivity:</b>	Human,Mouse,Rat

---

### Additional Information

<b>Observed MW:</b>	78kDa	<b>Conjugate:</b>	Unconjugated
<b>Calculated MW:</b>	65kDa	<b>Isotype:</b>	IgG

---

### Immunogen Information

**Background:** This gene encodes a component of the interleukin 1 receptor complex, which initiates signalling events that result in the activation of interleukin 1-responsive genes. Alternative splicing of this gene results in membrane-bound and soluble isoforms differing in their C-terminus. The ratio of soluble to membrane-bound forms increases during acute-phase induction or stress.

**Recommended Dilution:** WB,1:1000 - 1:5000

**Synonyms:** IL1R3; C3orf13; IL-1RAcP

**Purification Method:** Affinity purification

**Immunogen:** Recombinant fusion protein containing a sequence corresponding to amino acids 391-470 of human IL-1RAcP (NP\_002173.1).

**Storage:** Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.05% proclin300,0.05% BSA,50% glycerol,pH7.3.