

CAB22294

Product Information

Product SKU:	CAB22294	Gene ID:	-	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Arabidopsis thaliana, Zea mays, Triticum aestivum, Other (Wide Range Predicted)

Additional Information

Observed MW:	15kDa	Conjugate:	Unconjugated
Calculated MW:	14kDa	Isotype:	IgG

Immunogen Information

Background:

Recommended Dilution: WB, 1:500 - 1:1000

Synonyms: -

Purification Method: Affinity purification

Immunogen: A synthetic peptide corresponding to a sequence of Triticum aestivum ATG8.

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300, 50% glycerol, pH7.3.