

CAB22369

Product Information

Product SKU:	CAB22369	Gene ID:	8809	Size:	20uL, 100uL
Clone No:	ARC52215	Host Species:	Rabbit	Reactivity:	Human,Mouse

Additional Information

Observed MW:	80-105kDa	Conjugate:	Unconjugated
Calculated MW:	62kDa	Isotype:	IgG

Immunogen Information

Background: The protein encoded by this gene is a cytokine receptor that belongs to the interleukin 1 receptor family. This receptor specifically binds interleukin 18 (IL18), and is essential for IL18 mediated signal transduction. IFN-alpha and IL12 are reported to induce the expression of this receptor in NK and T cells. This gene along with four other members of the interleukin 1 receptor family, including IL1R2, IL1R1, ILRL2 (IL-1Rrp2), and IL1RL1 (T1/ST2), form a gene cluster on chromosome 2q. Alternatively spliced transcript variants encoding different isoforms have been found for this gene.

Recommended Dilution: WB,1:500 - 1:1000

Synonyms: CD218a; IL18RA; IL1RRP; CDw218a; IL-1Rrp; IL-18Ralpha; IL18Ralpha2; IL-18R-alpha; IL18R1/IL-18 R alpha

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 141-320 of human IL18R1/IL-18 R alpha (NP_003846.1).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.05% proclin300,0.05% BSA,50% glycerol,pH7.3.