

CAB22404

Product Information

Product SKU:	CAB22404	Gene ID:	-	Size:	20uL, 100uL
Clone No:	ARC54577	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat,Other (Wide Range Predicted)

Additional Information

Observed MW:	Refer to figures	Conjugate:	-
Calculated MW:	-	Isotype:	IgG

Immunogen Information

Background:

Recommended Dilution: DB,1:500 - 1:1000

Synonyms: m5C

Purification Method: Affinity purification

Immunogen: Chemical compounds corresponding to m5C.

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.05% proclin300,0.05% BSA,50% glycerol,pH7.3.