

CAB22411

Product Information

Product SKU:	CAB22411	Gene ID:	CAS:1867-73-8	Size:	20uL, 100uL
Clone No:	ARC5003-03	Host Species:	Rabbit	Reactivity:	Species independent

Additional Information

Observed MW:	Refer to figures	Conjugate:	-
Calculated MW:	-	Isotype:	IgG

Immunogen Information

Background:	Discovered in the 1970s, m6A is the most prevalent internal modification in polyadenylated mRNAs and long non-coding RNAs (lncRNAs) in higher eukaryotes. m6A is widely conserved among eukaryotic species that range from yeast, plants, flies to mammals, as well as among viral RNAs with a nuclear phase. The m6A-based modification is associated with a well-defined RNA motif, RRACH (R: A/G, H: A/C/U). As a representative of the epitranscriptome, m6A mRNA modifications participate in many vital activities in the cell, including stem cell self-renewal and differentiation, mRNA transcription, alternative splicing, nuclear export, translation, degradation, and microRNA processing. These processes determine the expression or inactivation of specific genes, which is vital for growth and development.(PMID: 30416848; PMID: 24662220; PMID: 30429466)
Recommended Dilution:	DB,1:500 - 1:2000 IF/ICC,1:50 - 1:200 meRIP,1:50 - 1:200
Synonyms:	N6-methyladenosine; m6A; N6-methyladenosine / m6A
Purification Method:	Affinity purification
Immunogen:	Chemical compounds corresponding to N6-methyladenosine / m6A.
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,0.05% BSA,50% glycerol,pH7.3.

Contact Details | Dublin, Ireland

Email: techsupport@assaygenie.com | **Web:** www.assaygenie.com

Copyright © 2020 Reagent Genie, All Rights Reserved. All information / detail is correct at time of going to print.