

## CAB22471

---

### Product Information

<b>Product SKU:</b>	CAB22471	<b>Gene ID:</b>	7048	<b>Size:</b>	20uL, 100uL
<b>Clone No:</b>	ARC57218	<b>Host Species:</b>	Rabbit	<b>Reactivity:</b>	Human,Mouse

---

### Additional Information

<b>Observed MW:</b>	Refer to figures	<b>Conjugate:</b>	Unconjugated
<b>Calculated MW:</b>	65kDa	<b>Isotype:</b>	IgG

---

### Immunogen Information

<b>Background:</b>	The protein encoded by this gene is a transmembrane protein that has a protein kinase domain, forms a heterodimeric complex with TGF-beta receptor type-1, and binds TGF-beta. This receptor/ligand complex phosphorylates proteins, which then enter the nucleus and regulate the transcription of genes related to cell proliferation, cell cycle arrest, wound healing, immunosuppression, and tumorigenesis. Mutations in this gene have been associated with Marfan Syndrome, Loeys-Deitz Aortic Aneurysm Syndrome, and the development of various types of tumors. Alternatively spliced transcript variants encoding different isoforms have been characterized.
<b>Recommended Dilution:</b>	IHC-P,1:50 - 1:200
<b>Synonyms:</b>	AAT3; FAA3; LDS2; MFS2; RIIC; LDS1B; LDS2B; TAAD2; TBR1; TBR-ii; TGFR-2; tbetaR-II; TGFbeta-R11; TGF beta Receptor II (TGFB2)
<b>Purification Method:</b>	Affinity purification
<b>Immunogen:</b>	A synthetic peptide corresponding to a sequence within amino acids 200-300 of human TGF beta Receptor II (TGFB2) (NP_003233.4).
<b>Storage:</b>	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.05% proclin300,0.05% BSA,50% glycerol,pH7.3.