

## CAB22549

---

### Product Information

<b>Product SKU:</b>	CAB22549	<b>Gene ID:</b>	22834	<b>Size:</b>	20uL, 100uL
<b>Clone No:</b>	-	<b>Host Species:</b>	Rabbit	<b>Reactivity:</b>	Human,Mouse,Rat

---

### Additional Information

<b>Observed MW:</b>	72kDa	<b>Conjugate:</b>	-
<b>Calculated MW:</b>	70kDa	<b>Isotype:</b>	IgG

---

### Immunogen Information

**Background:** Predicted to enable DNA-binding transcription factor activity, RNA polymerase II-specific and RNA polymerase II cis-regulatory region sequence-specific DNA binding activity. Predicted to be involved in regulation of transcription by RNA polymerase II. Predicted to be active in nucleus.

**Recommended Dilution:** WB,1:500 - 1:1000

**Synonyms:** -

**Purification Method:** Affinity purification

**Immunogen:** Recombinant fusion protein containing a sequence corresponding to amino acids 1-74 of human ZNF652 (NP\_055712.1).

**Storage:** Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.