

## CAB22550

---

### Product Information

<b>Product SKU:</b>	CAB22550	<b>Gene ID:</b>	23025	<b>Size:</b>	20uL, 100uL
<b>Clone No:</b>	-	<b>Host Species:</b>	Rabbit	<b>Reactivity:</b>	Human,Mouse,Rat

---

### Additional Information

<b>Observed MW:</b>	220kDa	<b>Conjugate:</b>	-
<b>Calculated MW:</b>	193kDa	<b>Isotype:</b>	IgG

---

### Immunogen Information

**Background:** This gene encodes a member of the UNC13 family. UNC13 proteins bind to phorbol esters and diacylglycerol and play important roles in neurotransmitter release at synapses. Single nucleotide polymorphisms in this gene may be associated with sporadic amyotrophic lateral sclerosis.

**Recommended Dilution:** WB,1:1000 - 1:5000

**Synonyms:** Munc13-1; Munc 13-1

**Purification Method:** Affinity purification

**Immunogen:** Recombinant protein of human Munc 13-1.

**Storage:** Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.