

CAB22551

Product Information

Product SKU:	CAB22551	Gene ID:	55288	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	75kDa	Conjugate:	Unconjugated
Calculated MW:	71kDa	Isotype:	IgG

Immunogen Information

Background: Predicted to enable GTP binding activity and GTPase activity. Involved in cellular homeostasis; mitochondrial outer membrane permeabilization; and mitochondrion transport along microtubule. Is integral component of mitochondrial outer membrane.

Recommended Dilution: WB,1:100 - 1:500 IF/ICC,1:50 - 1:200

Synonyms: ARHT1; MIRO1; MIRO-1; [KD Validated] RHOT1/MIRO1

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 310-590 of human RHOT1/MIRO1. (NP_060777.3).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.