

## CAB22554

---

### Product Information

<b>Product SKU:</b> CAB22554	<b>Gene ID:</b> 57535	<b>Size:</b> 20uL, 100uL
<b>Clone No:</b> -	<b>Host Species:</b> Rabbit	<b>Reactivity:</b> Human,Mouse,Rat

---

### Additional Information

<b>Observed MW:</b> 140kDa/	<b>Conjugate:</b> Unconjugated
<b>Calculated MW:</b> 111kDa	<b>Isotype:</b> IgG

---

### Immunogen Information

**Background:** Expression of this gene is induced by estrogen and the encoded protein has been characterized as a transmembrane protein. The encoded protein has been found in to correlate with survival in certain carcinomas (PMID: 21102415) and may be important for cellular response to stress (PMID: 21072319). Alternative splicing results in multiple transcript variants.

**Recommended Dilution:** WB,1:500 - 1:1000

**Synonyms:** EIG121; KIAA1324; ELAPOR1

**Purification Method:** Affinity purification

**Immunogen:** Recombinant protein of human ELAPOR1.

**Storage:** Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.