

CAB22568

Product Information

Product SKU:	CAB22568	Gene ID:	4738	Size:	20uL, 100uL
Clone No:	ARC58460	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	9kDa	Conjugate:	Unconjugated
Calculated MW:	9kDa	Isotype:	IgG

Immunogen Information

Background: Enables ubiquitin protein ligase binding activity. Acts upstream of or within protein neddylation. Located in cytosol and nucleoplasm. Biomarker of Parkinson's disease and malignant astrocytoma.

Recommended Dilution: WB,1:2000 - 1:20000

Synonyms: NEDD-8; NEDD8

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-81 of human NEDD8 (NP_006147.1).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.05% proclin300,0.05% BSA,50% glycerol,pH7.3.