

CAB22597

Product Information

Product SKU:	CAB22597	Gene ID:	27063	Size:	20uL, 100uL
Clone No:	ARC55780_ARC55783	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	36kDa	Conjugate:	Unconjugated
Calculated MW:	36kDa	Isotype:	IgG

Immunogen Information

Background: The protein encoded by this gene is localized to the nucleus of endothelial cells and is induced by IL-1 and TNF-alpha stimulation. Studies in rat cardiomyocytes suggest that this gene functions as a transcription factor. Interactions between this protein and the sarcomeric proteins myopalladin and titin suggest that it may also be involved in the myofibrillar stretch-sensor system.

Recommended Dilution: WB,1:500 - 1:1000

Synonyms: ALRP; CARP; C-193; CVARP; MCARP; bA320F15.2

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-130 of human ANKRD1 Rabbit PolymAb® (NP_055206.2).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.05% proclin300,0.05% BSA,50% glycerol,pH7.3.