

CAB22617

Product Information

Product SKU:	CAB22617	Gene ID:	3702	Size:	20uL, 100uL
Clone No:	ARC58362	Host Species:	Rabbit	Reactivity:	Human

Additional Information

Observed MW:	75kDa	Conjugate:	Unconjugated
Calculated MW:	72kDa	Isotype:	IgG

Immunogen Information

Background: This gene encodes an intracellular tyrosine kinase expressed in T-cells. The protein contains both SH2 and SH3 domains which are often found in intracellular kinases. It is thought to play a role in T-cell proliferation and differentiation.

Recommended Dilution: WB, 1:2000 - 1:6000

Synonyms: EMT; LYK; LPFS1; PSCTK2; ITK/EMT

Purification Method: Affinity purification

Immunogen: Recombinant protein of human ITK/EMT.

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300, 0.05% BSA, 50% glycerol, pH7.3.