

CAB22622

Product Information

Product SKU:	CAB22622	Gene ID:	8100	Size:	20uL, 100uL
Clone No:	ARC57206	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	95kDa	Conjugate:	Unconjugated
Calculated MW:	93kDa	Isotype:	IgG

Immunogen Information

Background: This gene encodes a member of the tetratrico peptide repeat (TPR) family. The encoded protein is involved in cilium biogenesis. Mutations of a similar gene in mouse can cause polycystic kidney disease. Several transcript variants encoding distinct isoforms have been identified for this gene.

Recommended Dilution: WB,1:500 - 1:1000

Synonyms: DAF19; TG737; TTC10; hTg737; D13S1056E; IFT88

Purification Method: Affinity purification

Immunogen: Recombinant protein of human IFT88.

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.05% proclin300,0.05% BSA,50% glycerol,pH7.3.