

CAB22637

Product Information

Product SKU:	CAB22637	Gene ID:	3344	Size:	20uL, 100uL
Clone No:	ARC58450	Host Species:	Rabbit	Reactivity:	Human

Additional Information

Observed MW:	Refer to figures	Conjugate:	Unconjugated
Calculated MW:	47kDa	Isotype:	IgG

Immunogen Information

Background: This gene encodes a forkhead domain binding protein and may function in the transcriptional regulation of the human T-cell leukemia virus long terminal repeat.

Recommended Dilution: WB, 1:500 - 1:1000

Synonyms: HTLF; N2

Purification Method: Affinity purification

Immunogen: Recombinant protein of human FOXN2.

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300, 0.05% BSA, 50% glycerol, pH7.3.