

## CAB22651

---

### Product Information

<b>Product SKU:</b>	CAB22651	<b>Gene ID:</b>	-	<b>Size:</b>	20uL, 100uL
<b>Clone No:</b>	ARC58400	<b>Host Species:</b>	Rabbit	<b>Reactivity:</b>	Human

---

### Additional Information

<b>Observed MW:</b>	60kDa	<b>Conjugate:</b>	Unconjugated
<b>Calculated MW:</b>	51kDa	<b>Isotype:</b>	IgG

---

### Immunogen Information

**Background:**

**Recommended Dilution:** WB, 1:2000 - 1:20000

**Synonyms:** -

**Purification Method:** Affinity purification

**Immunogen:** Recombinant protein of Enterovirus EV71 3D.

**Storage:** Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300, 0.05% BSA, 50% glycerol, pH7.3.