

CAB22655

Product Information

Product SKU:	CAB22655	Gene ID:	10841	Size:	20uL, 100uL
Clone No:	ARC58234	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	59kDa	Conjugate:	Unconjugated
Calculated MW:	59kDa	Isotype:	IgG

Immunogen Information

Background: The protein encoded by this gene is a bifunctional enzyme that channels 1-carbon units from formiminoglutamate, a metabolite of the histidine degradation pathway, to the folate pool. Mutations in this gene are associated with glutamate formiminotransferase deficiency. Alternatively spliced transcript variants have been found for this gene.

Recommended Dilution: WB,1:2000 - 1:9000

Synonyms: LCHC1; FTCD

Purification Method: Affinity purification

Immunogen: Recombinant protein of human FTCD.

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.05% proclin300,0.05% BSA,50% glycerol,pH7.3.