

CAB22656

Product Information

Product SKU:	CAB22656	Gene ID:	139341	Size:	20uL, 100uL
Clone No:	ARC57663	Host Species:	Rabbit	Reactivity:	Rat

Additional Information

Observed MW:	Refer to figures	Conjugate:	Unconjugated
Calculated MW:	17kDa	Isotype:	IgG

Immunogen Information

Background: This gene encodes a protein with a FUN14 superfamily domain. The function of the encoded protein is not known.

Recommended Dilution: IHC-P,1:100 - 1:500

Synonyms: FUNDC1

Purification Method: Affinity purification

Immunogen: A synthesized peptide derived from human FUNDC1.

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300,0.05% BSA,,50% glycerol,pH7.3.