

CAB22660

Product Information

Product SKU:	CAB22660	Gene ID:	6441	Size:	20uL, 100uL
Clone No:	ARC57711	Host Species:	Rabbit	Reactivity:	Mouse,Rat

Additional Information

Observed MW:	50kDa	Conjugate:	Unconjugated
Calculated MW:	38kDa	Isotype:	IgG

Immunogen Information

Background: The protein encoded by this gene is part of the innate immune response, protecting the lungs against inhaled microorganisms and chemicals. The encoded protein may also be involved in surfactant metabolism.

Recommended Dilution: WB,1:2000 - 1:20000

Synonyms: SP-D; PSP-D; SFTP4; COLEC7; SFTPD

Purification Method: Affinity purification

Immunogen: A synthesized peptide derived from human Surfactant protein D/SP-D.

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.05% proclin300,0.05% BSA,50% glycerol,pH7.3.