

CAB22669

Product Information

Product SKU:	CAB22669	Gene ID:	39105218	Size:	20uL, 100uL
Clone No:	ARC57815	Host Species:	Rabbit	Reactivity:	Human

Additional Information

Observed MW:	110kDa	Conjugate:	Unconjugated
Calculated MW:	150kDa	Isotype:	IgG

Immunogen Information

Background: S1 attaches the virion to the cell membrane by interacting with sialic acid-containing cell receptors, initiating the infection.; [Spike protein S1]: attaches the virion to the cell membrane by interacting with host receptor, initiating the infection.

Recommended Dilution: WB, 1:2000 - 1:9000

Synonyms: -

Purification Method: Affinity purification

Immunogen: A synthesized peptide derived from HCoV-OC43 Spike S1.

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300, 0.05% BSA, 50% glycerol, pH7.3.