

CAB22679

Product Information

Product SKU:	CAB22679	Gene ID:	27036	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human

Additional Information

Observed MW:	70kDa	Conjugate:	Unconjugated
Calculated MW:	51kDa	Isotype:	IgG

Immunogen Information

Background: Predicted to enable sialic acid binding activity. Predicted to be involved in cell adhesion. Predicted to be integral component of plasma membrane. Predicted to be active in plasma membrane.

Recommended Dilution: WB, 1:500 - 1:1000

Synonyms: p75; QA79; AIRM1; CD328; AIRM-1; CDw328; D-siglec; SIGLEC-7; SIGLECP2; SIGLEC19P; p75/AIRM1

Purification Method: Affinity purification

Immunogen: Recombinant protein of human CD328.

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.