

CAB22698

Product Information

Product SKU:	CAB22698	Gene ID:	54973	Size:	20uL, 100uL
Clone No:	ARC58804	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	68kDa	Conjugate:	Unconjugated
Calculated MW:	68kDa	Isotype:	IgG

Immunogen Information

Background: The Integrator complex contains at least 12 subunits and associates with the C-terminal domain of RNA polymerase II large subunit (POLR2A; MIM 180660) and mediates the 3-prime end processing of small nuclear RNAs U1 (RNU1; MIM 180680) and U2 (RNU2; MIM 180690). INTS11, or CPSF3L, is the catalytic subunit of the Integrator complex (Baillat et al., 2005 [PubMed 16239144]).

Recommended Dilution: WB, 1:2000 - 1:7000

Synonyms: RC68; INT11; RC-68; CPSF3L; CPSF73L; INTS11

Purification Method: Affinity purification

Immunogen: Recombinant protein of human INTS11.

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300, 0.05% BSA, 50% glycerol, pH7.3.