

CAB22739

---

## Product Information

<b>Product SKU:</b>	CAB22739	<b>Gene ID:</b>	39105218	<b>Size:</b>	20uL, 100uL
<b>Clone No:</b>	ARC57814	<b>Host Species:</b>	Rabbit	<b>Reactivity:</b>	Human

---

## Additional Information

<b>Observed MW:</b>	110kDa	<b>Conjugate:</b>	Unconjugated
<b>Calculated MW:</b>	150kDa	<b>Isotype:</b>	IgG

---

## Immunogen Information

**Background:** S1 attaches the virion to the cell membrane by interacting with sialic acid-containing cell receptors, initiating the infection.; [Spike protein S1]: attaches the virion to the cell membrane by interacting with host receptor, initiating the infection.

**Recommended Dilution:** WB, 1:2000 - 1:9000

**Synonyms:** -

**Purification Method:** Affinity purification

**Immunogen:** A synthesized peptide derived from human HCoV-OC43 Spike S1.

**Storage:** Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300, 0.05% BSA, 50% glycerol, pH 7.3.