

CAB22753

Product Information

Product SKU:	CAB22753	Gene ID:	1088	Size:	20uL, 100uL
Clone No:	ARC57403	Host Species:	Rabbit	Reactivity:	Human

Additional Information

Observed MW:	Refer to figures	Conjugate:	Unconjugated
Calculated MW:	38kDa	Isotype:	IgG

Immunogen Information

Background: Enables protein heterodimerization activity. Involved in heterophilic cell-cell adhesion via plasma membrane cell adhesion molecules. Located in cell surface and extracellular space. Biomarker of severe acute respiratory syndrome.

Recommended Dilution: IHC-P,1:500 - 1:1000

Synonyms: CD67; CGM6; CD66b; NCA-95; CD66b/CEACAM8

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 35-319 of human CD66b/CEACAM8.

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300,0.05% BSA,50% glycerol,pH7.3.