

CAB22755

Product Information

Product SKU:	CAB22755	Gene ID:	65980	Size:	20uL, 100uL
Clone No:	ARC58081	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	80kDa	Conjugate:	Unconjugated
Calculated MW:	67kDa	Isotype:	IgG

Immunogen Information

Background: Enables lysine-acetylated histone binding activity. Predicted to be involved in regulation of transcription by RNA polymerase II. Located in nucleoplasm. Part of SWI/SNF complex.

Recommended Dilution: WB,1:1000 - 1:5000

Synonyms: PRO9856; SMARCI2; LAVS3040; BRD9

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 274-597 of human BRD9 (NP_076413.3).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.05% proclin300,0.05% BSA,50% glycerol,pH7.3.