

## CAB22789

---

### Product Information

<b>Product SKU:</b>	CAB22789	<b>Gene ID:</b>	841395	<b>Size:</b>	20uL, 100uL
<b>Clone No:</b>	-	<b>Host Species:</b>	Rabbit	<b>Reactivity:</b>	Arabidopsis thaliana

---

### Additional Information

<b>Observed MW:</b>	-	<b>Conjugate:</b>	Unconjugated
<b>Calculated MW:</b>	39kDa	<b>Isotype:</b>	IgG

---

### Immunogen Information

**Background:** Identified as a protein that binds to abscisic acid response elements. May mediate transcriptional regulation of ABA responses.

**Recommended Dilution:** WB, 1:500 - 1:1000

**Synonyms:** abscisic acid responsive element-binding factor 1; AtABF1; F14J22.7; F14J22\_7; ABF1

**Purification Method:** Affinity purification

**Immunogen:** Recombinant protein of Arabidopsis thaliana ABF1.

**Storage:** Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.