

## CAB22792

---

### Product Information

<b>Product SKU:</b> CAB22792	<b>Gene ID:</b> 55726	<b>Size:</b> 20uL, 100uL
<b>Clone No:</b> ARC55758	<b>Host Species:</b> Rabbit	<b>Reactivity:</b> Human

---

### Additional Information

<b>Observed MW:</b> 80kDa	<b>Conjugate:</b> Unconjugated
<b>Calculated MW:</b> 80kDa	<b>Isotype:</b> IgG

---

### Immunogen Information

**Background:** Involved in regulation of mitotic cell cycle. Acts upstream of or within centrosome localization; mitotic spindle organization; and protein localization to nuclear envelope. Located in cytoplasm and nuclear body.

**Recommended Dilution:** WB, 1:10000 - 1:80000

**Synonyms:** ASUN; GCT1; NET48; Mat89Bb; SPATA30; C12orf11; ASUN/INTS13

**Purification Method:** Affinity purification

**Immunogen:** Recombinant protein of human ASUN/INTS13.

**Storage:** Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300, 0.05% BSA, 50% glycerol, pH7.3.