

CAB22866

Product Information

Product SKU:	CAB22866	Gene ID:	3783729	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human

Additional Information

Observed MW:	100kDa	Conjugate:	Unconjugated
Calculated MW:	56kDa	Isotype:	IgG

Immunogen Information

Background: BRRF2 is a capsid protein, but it can be secreted into the bloodstream, giving it the opportunity to be a diagnostic or immunotherapy target. The Epstein-Barr virus BRRF2 gene product is involved in viral progeny production.

Recommended Dilution: WB, 1:500 - 1:1000

Synonyms: Tegument protein BRRF2; BRRF2

Purification Method: Affinity purification

Immunogen: Recombinant protein of Epstein-Barr virus (strain B95-8) (HHV-4) BRRF2.

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH 7.3.