

CAB22868

Product Information

Product SKU: CAB22868	Gene ID: 607	Size: 20uL, 100uL
Clone No: -	Host Species: Rabbit	Reactivity: Mouse

Additional Information

Observed MW: 149kDa	Conjugate: -
Calculated MW: 149kDa	Isotype: IgG

Immunogen Information

Background: BCL9 is associated with B-cell acute lymphoblastic leukemia. It may be a target of translocation in B-cell malignancies with abnormalities of 1q21. Its function is unknown. The overexpression of BCL9 may be of pathogenic significance in B-cell malignancies.

Recommended Dilution: WB, 1:500 - 1:1000

Synonyms: LGS; Bcl9

Purification Method: Affinity purification

Immunogen: Recombinant protein of human Bcl9.

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.