

CAB22873

Product Information

Product SKU:	CAB22873	Gene ID:	9422	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human

Additional Information

Observed MW:	71kDa	Conjugate:	Unconjugated
Calculated MW:	71kDa	Isotype:	IgG

Immunogen Information

Background: This gene encodes a zinc finger protein and belongs to the krueppel C2H2-type zinc-finger protein family. Zinc finger proteins are often localized in the nucleus, bind nucleic acids, and regulate transcription.

Recommended Dilution: WB, 1:500 - 1:2000

Synonyms: -

Purification Method: Affinity purification

Immunogen: Recombinant protein of human ZNF264.

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.