

CAB22991

Product Information

Product SKU:	CAB22991	Gene ID:	51673	Size:	20uL, 100uL
Clone No:	ARC58837	Host Species:	Rabbit	Reactivity:	Mouse,Rat

Additional Information

Observed MW:	19kDa	Conjugate:	Unconjugated
Calculated MW:	19kDa	Isotype:	IgG

Immunogen Information

Background: Enables tubulin binding activity. Involved in decidualization and microtubule bundle formation. Colocalizes with microtubule bundle and perinuclear region of cytoplasm.

Recommended Dilution: WB,1:500 - 1:1000

Synonyms: p20; CGI-38; TPPP/p20; p25gamma; TPPP3

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-176 of human TPPP3(NP_057048.2).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.05% proclin300,0.05% BSA,50% glycerol,pH7.3.